

My Memorable Mistakes

Sona Damechani

Hematopathology Lab

TMH –Parel.

Introduction..

- A laboratory error is defined as any defect that occurs during the entire testing process, from ordering tests to reporting results, that in any way influences the quality of laboratory services.
- Any error during the laboratory testing process can affect patient care, including delay in reporting, unnecessary redraws, misdiagnosis, and improper treatment.
- Sometimes, these errors may even be fatal.

Total testing process

```
graph TD; A[Total testing process] --> B[Analytical specimen testing.]; B --> C[Pre analytical]; B --> D[Post analytical]; C --> E[1. Test requisition]; C --> F[2. Collection]; C --> G[3. Transport]; C --> H[4. Accessioning]; D --> I[1. Reporting]; D --> J[2. Interpretation]; D --> K[3. Follow up]; D --> L[4. Storage]; D --> M[5. Retesting];
```

Analytical specimen testing.

Pre analytical

1. Test requisition
2. Collection
3. Transport
4. Accessioning

Post analytical

1. Reporting
2. Interpretation
3. Follow up
4. Storage
5. Retesting

Error-1

- Sample received from ICU for coagulation profile, test performed in the lab showed -----

RESULTS

Test	Result	Unit	Normal Range
Prothrombin Time	22	sec	12.7-15
MNPT	13	sec	
INR	1.76		0.8-1.2
APTT	60.8	sec	24 - 35.3
control	30	sec	

Action taken as per lab policy...

- Delta check of previous test results
- check clinical history

If abnormal Find out the cause

If previous results are normal Repeat assay

(In this case previous result is normal)

Repeat sample collected shows results as follows

TEST	RESULT	UNITS
PT	13.2	SECs
MNPT	13.5	SECs
APTT	31.0	SECs
CONTROL	29.5	SECs

This results are normal

- Possibilities . . .

1. Clinical conditions like patient suffering from bleeding, extensive transfusion, swelling, history of hematoma etc.
(Not suggestive of the abnormal results)

2. Patient on anticoagulant .
(NO)

3. c/o post operative.
(NO)

- **Pre analytical**


- Improper blood to anticoagulant ratio.
- Sample quality
- Central line collection

- **Analytical error**

- Reagent Deteriorated , mixing of reagents, partial aspiration of reagent or sample, problem with Analyzer.

- **Sample Quality**
- **Sample Quantity**
- **Reagents**
- **Instrument issues**
- **Clinical conditions of patient**

On enquiring about the collection timings and site it was made clear that it was **central line collection where, heparin is flushed through the line before collecting the sample.**



Proper procedure is to be followed while collecting sample from central line

Tips and warning

- Avoid blood collection from heparin-flushed lines to avoid heparin contamination and specimen dilution.
- If heparin line must be used, flush line with 5-mL saline and discard the first 5 mL of blood (or 6 line volumes) of blood.
- Underfilling or overfilling results in an imbalance in the blood to anticoagulant ratio and can produce artificially prolonged or shortened clotting times, respectively.

Error-2

- A 35 years /F Sample received for complete blood count -----

The screenshot displays a hematology analyzer software interface. The main window shows a table of test results with the following data:


Test	CRes	F	D	Prev.Run	F	D	Time	Prev.Res
WBC	4.40	**						
RBC	3.58							
HGB	10.0							
HCT	32.8							
MCV	91.6							
MCH	28.0							
MCHC	30.5							
CHCM	31.3							
CH	28.6							
RDW	20.4							
RDW	2.17							
PLT	16	**						
MPV	9.1	**						
%NEUT	51.0	**						
%LYMPH	12.2	**						
%MONO	29.8	**						
%EOS	0.1	**						
%BASO	0.3	**						
%LUC	6.6	**						
%NRBC	*****	**						
#NEUT	2.25	**						
#LYMPH	0.54	**						
#MONO	1.31	**						

The PLT (Platelet Count) value of 16 is circled in red. The morphology flags for HYPD (Hyperplatelets) and PLTCUM (Platelet Clumps) are also circled in red. The sample/system flags section shows the following data:

Sample/System flags	NRCELL	WBC	%NEUT	%LYMPH	%MONO	%EOS	%LUC	%NRBC
NRCELL	WBC	%NEUT	%LYMPH	%MONO	%EOS	%LUC	%NRBC	
PLT-CL	PLT	MPV	WBCP	PDW	PCT			
PX-NV	%NEUT	%LYMPH	%MONO	%EOS	%LUC			
WBC-CE	WBC	WBCu						

The software interface also includes a central window with flow cytometry plots and a right-hand panel with an 'Exit' button.

Peripheral Blood Smear – WRIGHT STAINING



PLATELET SEEN IN CLUMPLS

Causes of Platelet clumping in EDTA Sample

- Slow collection
- Improper mixing
- EDTA induced psuedothrombocytopenia.

Corrective action taken


Repeat collection was advised in

1) EDTA



2) Na-Citrate anticoagulant.

NOTE: Na-Citrate collection is only for Platelet count & not for reporting of other CBC parameters.

Citrate sample run on analyzer for Platelet count



EDTA VS CITRATE SAMPLE PERIPHERAL BLOOD SMEAR



Calculation for Plt count of Citrate sample

$$\text{Plt count} = (\text{citrate plt count} \times 1.1) \times 10^9/\text{L}$$

$$\text{Plt count} = (254 \times 1.1) \times 10^9/\text{L}$$

$$\text{Plt count} = 280 \times 10^9/\text{L}$$

Lombarts AJ, de Kieviet W. Recognition and prevention of pseudothrombocytopenia and concomitant pseudoleukocytosis. Am J Clin Pathol. 1988;89:634–639.

CASE 1 - EDTA induced Pseudothrombocytopenia

- ❖ EDTA- induced pseudothrombocytopenia is a common laboratory phenomenon which gives spurious low platelet count.
- ❖ It is often overlooked because blood smears are not evaluated visually in routine practice and histograms as well as warning flags of hematology analyzers are not interpreted correctly.
- ❖ Awareness of typical PLT and WBC patterns may alert to the presence of EDTA-PTCP in routine laboratory practice helps to avoid unnecessary investigations and over-treatment.

Conclusion

It is important to carry out our work in lab with utmost care and alertness.

- To follow a technically correct methodology from time of collection to reporting the results ,so that the results are precise, accurate and reliable.
- To avoid erroneous results due to technical errors.

Prevention is Better than cure!



'It costs less to prevent a problem than it does to correct it'

A formal quality system in the laboratory should prevent mistakes by means of:

- quality assurance measures
- quality control of the analytical results
- thorough documentation of the system
- efficient maintenance of records
- regular audits of all aspects of the system

**IT'S NOT A MISTAKE
TO MAKE A MISTAKE,
BUT IT'S A MISTAKE
TO REPEAT THE
SAME MISTAKE.**



Thank
you